Store MicroSureH₂O_s vials in a horizontal position and above freezing. If colour separation occurs in the ampoule, shake vial prior to testing. Dispose used glass test vials in a 'Sharps' container. Do not open or empty vials for disposal.

INSTRUCTIONS FOR USE
(One water test is required for each water line)

Step 1. Collect water sample - Remove all handpieces and attachments from water lines. Turn water volume on the dental unit to high. Flush water lines for a minimum of 2 minutes. Label a clean cup with operatory and water line location. Put on gloves and lay out tray paper to prepare testing area. Collect a minimum of 30 ml of sample water from a water line into a clean cup.

Step 2. Filter water sample - Attach the filter cartridge to the syringe. Pull up the plunger to the 10 ml mark on the syringe then push fully down on the plunger to the 0ml mark to ensure clear flow thru the filter. Place filter tip into the test water ensuring that the tip remains immersed at all times to avoid trapping air in the filter. Slowly draw in 10 ml of sample water. This may take 60-90 seconds to complete. If a 10 ml sample is collected, remove the filter from the syringe and empty the water from the syringe as it is no longer required. **NOTE: If it becomes difficult to collect the 10ml sample, the filter is saturated with bacteria and/or debris. Proceed to shock treatment before starting a new test.**

Step 3. Wash the sample - Reattach the filter cartridge to the syringe. Remove the blue lid from the vial with the clear solution. Position the tip of the filter cartridge into the vial, and draw all the solution into the syringe. Washing with distilled water/saline will ensure your sample is free of interfering chemicals and will yield accurate results.

Step 4. Perform the test - Remove the cap of the pink reagent vial and place the cap face up so that the inside is not in contact with any surface. Do not touch the inside of the vial or the inside of the cap as this could result in a contaminated sample. Slowly press the plunger down to add 5 drops of the sample to the vial with the pink reagent. Recap the vial, and shake to mix the two liquids.

Step 5. Label this test vial appropriately, this is important when testing multiple water lines. Let the test vial sit at room temperature for one hour, then observe and note colour.

NOTE: (Optional) Incubate - If the vial fits into your dry-block incubator, you can incubate the vial for 15-30 minutes to achieve faster results, then observe and note colour.

Step 6. Observe Colour using a white backdrop & Record Results (SEE COLOUR CHART ON BACK OF BOX)

Go to amrakhealthcare.com/resources-2 for additional tips and documents.

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